



VAX-ID® a reliable patient-friendly intradermal injection adapter

One solution fits all people and accommodates different medicinal products



Controlled depth




Easy and safe



Low in pain

Needle Gauge    Average Injection Depth

 <b>32G</b>	0.5mm
 <b>30G</b>	0.6mm
 <b>27G</b>	0.8mm

VAX-ID® is compatible with syringes like:

HSW, B.Braun	<b>(Non-prefilled)</b>
Stevanato, Symatese, Gerresheimer, ...	<b>(Prefilled)</b>

Built for the applications that depend on precision, safety and patient comfort



**THERAPEUTIC**

Cancer and Immuno therapies,  
Dermatology



**PREVENTIVE**

Vaccination strategies, e.g. Rabies



**DIAGNOSTIC**

Tuberculin Skin Test (TST),  
Lymphatic tracers

**PARTNER WITH IDEVAX!**



We welcome partnerships with

- Preclinical and Clinical Academic & Industry Research
- Biotech & Pharma
- Dermatology Clinics
- Travel vaccination, TB Centers & Health services
- Public & Global health agencies & NGOs
- Distributors

Learn more: [www.idevax.com](http://www.idevax.com)

